



Low Volume Sample Fixation for Single Cell RNA Sequencing

Protocol

Legal Notices

Document 1020807 Rev E, May 2025
© 2025 Scale Biosciences, Inc.

Scale Biosciences
5601 Oberlin Dr., Suite 110
San Diego, CA 92121, United States
www.scale.bio
support@scale.bio

Scale Biosciences, Inc (“ScaleBio”). All rights reserved. No part of this document may be reproduced, distributed, or transmitted in any form or by any means, electronic or mechanical, including photocopying, recording, or any information storage and retrieval system, without the prior written permission of ScaleBio. This document is provided for information purposes only and is subject to change or withdrawal by ScaleBio at any time.

Disclaimer of Warranty

TO THE EXTENT PERMITTED BY APPLICABLE LAW, SCALEBIO PROVIDES THIS DOCUMENT “AS IS” WITHOUT WARRANTY OF ANY KIND, INCLUDING WITHOUT LIMITATION, ANY IMPLIED WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR NONINFRINGEMENT. IN NO EVENT WILL SCALEBIO BE LIABLE TO YOU OR ANY THIRD PARTY FOR ANY LOSS OR DAMAGE, DIRECT OR INDIRECT, FROM THE USE OF THIS DOCUMENT, INCLUDING WITHOUT LIMITATION, LOST PROFITS, LOST INVESTMENT, BUSINESS INTERRUPTION, GOODWILL, OR LOST DATA, EVEN IF SCALEBIO IS EXPRESSLY ADVISED IN ADVANCE OF THE POSSIBILITY OF SUCH LOSS OR DAMAGE. Any warranties applicable to the ScaleBio products are set forth in the Terms and Conditions accompanying such product and such Terms and Conditions are not modified in any way by the terms of this notice.

Trademark Information

ScaleBio may make reference to products or services provided by other companies using their brand names or company names solely for the purpose of clarity, and does not assert any ownership rights over those third-party marks or names. Images were created with BioRender.com

Patent Information

ScaleBio products may be covered by one or more patents as indicated at <https://scale.bio/legal-notice/>.

Terms and Conditions

The use of the ScaleBio products described herein is subject to ScaleBio’s Terms and Conditions that accompany the product, or such other terms as have been agreed to in writing between ScaleBio and the user.

Intended Use

All products and services described herein are intended FOR RESEARCH USE ONLY and NOT FOR USE IN DIAGNOSTIC PROCEDURES.

Table of Contents

Required Materials	4
Best Practices	6
Assay Introduction	7
Cell Preparation	8
Cell Fixation	9
Revision History	12

Required Materials

Consumables and Reagents Manufactured by ScaleBio

ScaleBio Sample Fixation Kit (PN 2020001) Consumables and Reagents.

Kit Module	Consumable	Part Number	Quantity	Storage Temperature
Module A (PN 2020002)	Wash Buffer	202100001	12	-20°C
Module B (PN 2020003)	Fixation Reagent	202110001	12	4°C

Consumables and Reagents Manufactured by Other Vendors

Consumable or Reagent	Supplier	Part Number
1X PBS without calcium or magnesium	Various	Various
DEPC (purity $\geq 96\%$) ^a	Millipore Sigma	D5758-25ML
DMSO (anhydrous) ^a	Thermo Fisher	D12345
Methanol (purity $\geq 99.9\%$) ^a	Fisher Scientific	A412-500
Pipette tips (nuclease-free, filtered, low retention for P1000, P200, P20) ^b	Various	Various
Pipette tips (nuclease-free, filtered, wide bore for P1000) ^b	Various	Various
5-mL DNA LoBind tubes ^b	Eppendorf	0030108310
1.5-mL DNA LoBind tubes ^b	Eppendorf	0030108418
Cell strainer (optional) ^c	VWR	10032-802
Cell counting dye [before fixation: trypan blue, after fixation: AO/PI, YOYO-1]	Various	Various

a. Other vendors may be used if reagent formulation remains the same.

b. Required for best assay performance.

c. Type and filter size appropriate for sample (ie, cell diameter and cell suspension volume).

Recommended Equipment

Item	Supplier	Part Number
Centrifuge with temperature control (1.5-mL tubes)	Various	Various
Vortex mixer	Various	Various
Pipettes (P1000, P200, P20, P10)	Various	Various
Cell counter	Various	Various
Chemical fume hood	Various	Various

Best Practices

For general laboratory practices:

- Calibrate and service pipettes every 12 months to ensure accurate sample volume transfer at each step.
- Store all reagents at the storage conditions recommended by the supplier.
- Thaw all reagents on ice, unless otherwise specified.
- Unless otherwise specified, vortex reagents.
- Open Fixation Reagent packing in a chemical fume hood.
- Handle Fixation Reagent, DEPC, and methanol in a chemical fume hood.
- Never reuse pipette tips or tubes.
- Use wide-bore tips for pipetting cell suspensions.
- Keep pipette tip boxes, reagent containers, and sample tubes closed when not in use.
- Wear suitable protective clothing, eyewear, and gloves.

For RNase-free sample processing:

- Use low-retention, RNase-free pipette tips and low-binding reaction tubes to prevent adsorption to plastic surfaces.
- Routinely wipe work surfaces with RNase AWAY to remove RNases, and with a 10% bleach cleaning solution to remove DNA amplicon contaminants.
- Wear disposable gloves and change them frequently.

Assay Introduction

The ScaleBio™ Sample Fixation Kit (PN 2020001) is intended for fixing cells prior to use with QuantumScale Single Cell RNA Kits or the ScaleBio Single Cell RNA Sequencing Kit v1.1. This protocol allows fixation of samples containing between 100,000 up to 1 million cells (or nuclei). This protocol makes handling of low cell (or nuclei) numbers easier and reduces sample loss. We recommend performing pilot fixation(s) to evaluate the recovery of samples with low input using this protocol, to determine if yields in a prospective sample will be of sufficient input for Library Preparation. Fixation can be performed at different timepoints, as samples can be stored for up to 12 months at -80°C prior to use.

Cell Preparation

Review the following table to prepare reagents for starting this section.

Source	Material
Other vendors	<ul style="list-style-type: none"> • 1X PBS (without calcium or magnesium) • 1.5-mL DNA LoBind tubes • Pipette tips (nuclease-free, filtered, wide bore for P1000) • Cell counting dye • Cell counter • Cell strainer (optional)



IMPORTANT: Use wide-bore pipette tips, and gentle pipette-mixing, for preparing cell suspensions to maintain sample quality. Cells must be washed with 1X PBS and resuspended in a maximum volume of 50 μ L prior to fixation.

Before you Begin

- Place 1X PBS on ice.
- Start with sufficient cells to obtain 100,000 to 1 million cells after washing and counting, assuming cell loss throughout the process.

Procedure

1. Obtain cells from culture, tissue dissociation, or thaw if frozen, and place **on ice**.
2. Transfer the cell suspension into a 1.5-mL DNA LoBind tube.
3. Centrifuge the cells at **500 x g** for **5 minutes** at **4°C**.



NOTE: Nuclei or small cells can be centrifuged for 8–10 minutes to maximize recovery.

4. Carefully remove the supernatant without disturbing the pellet.
5. Add **500 μ L** ice cold 1X PBS and resuspend the cells by flicking.
6. Centrifuge the cells at **500 x g** for **5 minutes** at **4°C**.
7. Carefully remove the supernatant without disturbing the pellet.
8. Resuspend cells with **50 μ L** of 1X PBS and place on ice.
9. Determine the concentration of the cell suspension using cell counting equipment. For accurate cell counting, use ≥ 2 μ L of cell suspensions and appropriate dilution factors recommended for your cell counting method.
10. (Optional) If cell suspensions appear clumpy, strain cells to a new 1.5-mL tube through a cell strainer and repeat [step 9](#) to recount the cell suspension. **Pipetting may be used in place of cell straining to reduce cell loss.**

Cell Fixation

Review the following table to prepare reagents before starting this section.


Source	Material	Take From	Place At	Brief Vortex	Brief Spin
Module A	Wash Buffer	-20°C	On ice	✗	✗
Module B	Fixation Reagent	4°C	RT	✗	✓
Other Vendors	1X PBS	—	On ice	—	—
	100% Methanol	—	On ice	—	—
	DMSO	—	RT	—	—
	DEPC	—	On ice	—	—
	Cell counting dye				
	1.5-mL DNA LoBind tubes				
	5-mL DNA LoBind tubes				
	Chemical fume hood				

Before you Begin

- Place single-cell suspensions with 100,000 to 1 million cells in 50 µL of 1X PBS in a 1.5-mL DNA LoBind tube on ice.
- Once thawed, invert the Wash Buffer to ensure it is fully mixed, then place on ice.
- Bring centrifuges (for 1.5-mL or 5-mL tubes) to 4°C.

Procedure

- Once the Fixation Reagent tube has been equilibrated to **room temperature**, briefly spin down. Fixation Reagent is lyophilized at the bottom of the tube and appears as a white pellet.



white pellet

- At **room temperature** and in a chemical fume hood, add **50 µL** of DMSO to the Fixation Reagent tube.
- Vortex at high speed with intermittent brief spins to dissolve the Fixation Reagent. This may take up to several minutes. Ensure that the pellet is fully dissolved before proceeding.
- Briefly spin down the Fixation Reagent tube.

5. **On ice**, prepare Complete Fixation Solution by combining in a 5-mL DNA LoBind tube the reagents in the specified order according to [Table 1](#). This is enough for the fixation of 1–6 samples.

Table 1: Complete Fixation Solution for 1–6 Samples

	Reagent	Volume (μL)
Step 1	Ice cold 100% methanol	2000
	Reconstituted Fixation Reagent	50
Vortex for 10 seconds. Use within 6 hours. Immediately before use, add the following:		
Step 2	DEPC	20
	Total volume	2070

6. Briefly vortex to mix.
7. Add **200 μL** of the Complete Fixation Solution **slowly** to the 1.5-mL DNA LoBind tube containing cells.
8. Gently pipette up and down three times to mix.



CAUTION: Rapid addition of Complete Fixation Solution or forceful pipette mixing may negatively impact fixation efficiency and data quality. Ensure that proper care is taken in pipetting steps.

9. Incubate **on ice** for **15 minutes**.
10. Add **500 μL** of Wash Buffer **slowly** to the tube and gently flick the tube to mix.
11. Centrifuge the tube at **500 x g** for **5 minutes** at **4°C**.



NOTE: Nuclei or small cells can be centrifuged for 8–10 minutes to maximize recovery.

12. Carefully remove the supernatant without disturbing the pellet, leaving ~100 μL of residual volume.
13. Add **500 μL** of Wash Buffer and resuspend cells by flicking the tube.
14. Centrifuge the tube at **500 x g** for **5 minutes** at **4°C**.
15. Carefully remove the supernatant without disturbing the pellet, leaving **~20–30 μL** of residual volume, as shown.



NOTE: If the fixation protocol started with 100,000–200,000 cells, skip [step 16](#) to avoid excessive dilution of the sample.

16. Add Wash Buffer for a total combined volume of **~50 μL**.

17. Determine the concentration of the cell suspension using cell counting equipment. For accurate cell counting, use ≥ 2 μL of cell suspensions and appropriate dilution factors recommended for your cell counting method.
18. Store the fixed cells, or proceed to QuantumScale Single Cell RNA or ScaleBio Single Cell RNA Sequencing Library Preparation.



Safe stopping point. The fixed cell suspension can be stored at -80°C for up to 12 months before proceeding with Library Preparation.

Revision History

Document, Revision	Revisions Date	Description of Change
1020807 Rev E	May 2025	Updated for QuantumScale
1020807 Rev D	March 2024	Initial release