



Technical Note

Vector Compatibility with CRISPR Kit

For Research Use Only.

Legal Notices

Document 1097102, Rev A, Jul 2024
© 2024 Scale Biosciences, Inc.

3210 Merryfield Row
San Diego, CA 92121, United States
<https://scale.bio/>
support@scale.bio

Scale Biosciences, Inc (“ScaleBio”). All rights reserved. No part of this document may be reproduced, distributed, or transmitted in any form or by any means, electronic or mechanical, including photocopying, recording, or any information storage and retrieval system, without the prior written permission of ScaleBio. This document is provided for information purposes only and is subject to change or withdrawal by ScaleBio at any time.

Disclaimer of Warranty:

TO THE EXTENT PERMITTED BY APPLICABLE LAW, SCALEBIO PROVIDES THIS DOCUMENT “AS IS” WITHOUT WARRANTY OF ANY KIND, INCLUDING WITHOUT LIMITATION, ANY IMPLIED WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR NONINFRINGEMENT. IN NO EVENT WILL SCALEBIO BE LIABLE TO YOU OR ANY THIRD PARTY FOR ANY LOSS OR DAMAGE, DIRECT OR INDIRECT, FROM THE USE OF THIS DOCUMENT, INCLUDING WITHOUT LIMITATION, LOST PROFITS, LOST INVESTMENT, BUSINESS INTERRUPTION, GOODWILL, OR LOST DATA, EVEN IF SCALEBIO IS EXPRESSLY ADVISED IN ADVANCE OF THE POSSIBILITY OF SUCH LOSS OR DAMAGE. Any warranties applicable to the ScaleBio products are set forth in the Terms and Conditions accompanying such product and such Terms and Conditions are not modified in any way by the terms of this notice.

Trademark Information:

ScaleBio may make reference to products or services provided by other companies using their brand names or company names solely for the purpose of clarity, and does not assert any ownership rights over those third-party marks or names. Images were created with BioRender.com

Patent Information:

ScaleBio products may be covered by one or more patents as indicated at: <https://scale.bio/legal-notice/>

Terms and Conditions:

The use of the ScaleBio products described herein is subject to ScaleBio's Terms and Conditions that accompany the product, or such other terms as have been agreed to in writing between ScaleBio and the user.

Intended Use:

All products and services described herein are intended FOR RESEARCH USE ONLY and NOT FOR USE IN DIAGNOSTIC PROCEDURES.

Introduction

The ScaleBio™ CRISPR Guide Enrichment Kit enables the capture of CRISPR guide sequences from CROP-seq style vectors through specific RT primers contained on the CRISPR RT Barcode Plate. An enrichment step after final cleanup of the RNA library uses a primer against the human U6 promoter to amplify the CRISPR guides. This Technical Note details the oligo sequences used in both the Cas9 CRISPR RT Barcode Plate and the CRISPR Enrichment PCR Module and outlines considerations taken to ensure vector compatibility.


Vector Compatibility

CROP-seq integrates CRISPR/Cas9 genome editing with single cell RNA sequencing to enable high-throughput functional genomics screening at the single cell level (1). This method employs vectors that embed a copy of the CRISPR guide target sequence within a larger mRNA sequence (sgRNA) downstream of a known promoter. This allows for the capture of the guide sequence alongside that of the RNA library. To confirm compatibility of your plasmid with the ScaleBio CRISPR Guide Enrichment Kit please verify your vector contains the following:

1. Guide RNA must be embedded in a longer transcript (see “Example Vector”)
2. Guide RNA scaffold should contain the RT capture primer sequence* (see “CRISPR Kit Oligos – Cas9-specific RT primer”)
3. Guide RNA must be downstream of the human U6 promoter* (see “CRISPR Kit Oligos – CRISPR Enrichment Module”)

*Note that if the human U6 sequence or scaffold primer sequence is not found in the vector/target region, custom primers can be designed to enable capture of other sequences. Please contact support@scale.bio for more information.

Example Vector



Compatible Vector Examples:

Vector	Addgene Link
CROPseq-Guide-Puro	86708
CROPseq-puro-v2	127458
CROP-seq-opti	106280

CRISPR Kit Oligos

Cas9-specific CRISPR RT Primer:

RT Primer structure:

[5' Phosphate Overhang] - [UMI] - [RT Barcode] - [cas9 Scaffold Capture*]

Example sequence for ordering oligos:

/5Phos/CAGAGC NNNNNNNN RT Barcode ACTTTTCAAGTTGATAACGGACTAGCCTTATTT

*Customization:

Note that the “cas9 Scaffold Capture” sequence can be replaced with any other sequence to prime other alternative Cas enzyme scaffolds, barcodes, or transcript regions. Please reach out to support@scale.bio for more information.

CRISPR Enrichment Module (downstream PCR primer):

U6-specific PCR Primer Sequence* :

P7 - i7 Index - Read2 - CTTGTGGAAGGACGAAACACC

*Customization:

Note that this primer can be replaced with another sequence if U6 is not found in the plasmid or region of interest. Please reach out to support@scale.bio for more information.

References

1. Datlinger P, Rendeiro AF, Schmidl C, Krausgruber T, Traxler P, Klughammer J, Schuster LC, Kuchler A, Alpar D, Bock C. Pooled CRISPR screening with single-cell transcriptome readout. Nat Methods. 2017 Mar;14(3):297-301.

Document Revision History

Revision	Revision Date	Document ID	Changes
Rev A	Jul 2024	1097102	Initial release.