

Technical Note

Guidelines for Pooling QuantumScale and 10x Genomics Libraries

For Research Use Only.

Legal Notices

Document 1352112, Rev A, Jun 2025

© 2025 Scale Biosciences, Inc.

5601 Oberlin Dr., Suite 110

San Diego, CA 92121

<https://scale.bio/>

support@scale.bio

Scale Biosciences, Inc (“ScaleBio”). All rights reserved. No part of this document may be reproduced, distributed, or transmitted in any form or by any means, electronic or mechanical, including photocopying, recording, or any information storage and retrieval system, without the prior written permission of ScaleBio. This document is provided for information purposes only and is subject to change or withdrawal by ScaleBio at any time.

Disclaimer of Warranty:

TO THE EXTENT PERMITTED BY APPLICABLE LAW, SCALEBIO PROVIDES THIS DOCUMENT “AS IS” WITHOUT WARRANTY OF ANY KIND, INCLUDING WITHOUT LIMITATION, ANY IMPLIED WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR NONINFRINGEMENT. IN NO EVENT WILL SCALEBIO BE LIABLE TO YOU OR ANY THIRD PARTY FOR ANY LOSS OR DAMAGE, DIRECT OR INDIRECT, FROM THE USE OF THIS DOCUMENT, INCLUDING WITHOUT LIMITATION, LOST PROFITS, LOST INVESTMENT, BUSINESS INTERRUPTION, GOODWILL, OR LOST DATA, EVEN IF SCALEBIO IS EXPRESSLY ADVISED IN ADVANCE OF THE POSSIBILITY OF SUCH LOSS OR DAMAGE. Any warranties applicable to the ScaleBio products are set forth in the Terms and Conditions accompanying such product and such Terms and Conditions are not modified in any way by the terms of this notice.

Trademark Information:

ScaleBio may make reference to products or services provided by other companies using their brand names or company names solely for the purpose of clarity, and does not assert any ownership rights over those third-party marks or names. Images were created with BioRender.com

Patent Information:

ScaleBio products may be covered by one or more patents as indicated at:

<https://scale.bio/legal-notice/>

Terms and Conditions:

The use of the ScaleBio products described herein is subject to ScaleBio's Terms and Conditions that accompany the product, or such other terms as have been agreed to in writing between ScaleBio and the user.

Intended Use:

All products and services described herein are intended FOR RESEARCH USE ONLY and NOT FOR USE IN DIAGNOSTIC PROCEDURES.

Introduction

Pooling Quantum libraries with 10x Genomics libraries on a single sequencing run is technically possible, but requires careful planning and comes with important caveats. The following guidance is based on an internal assessment by the Scale Bio technical team and is provided for informational purposes only. Performance of this approach has not been fully validated in production settings, and users should proceed cautiously and in close coordination with their bioinformatics and sequencing teams.

Disclaimer: This guidance is **hypothetical** and reflects a feasibility assessment only. It is not an officially supported configuration, and Scale Bio cannot guarantee sequencing or downstream analysis performance when Quantum and 10x libraries are pooled.

Pooling Quantum Libraries with 10x Genomics Libraries

Read Configuration Considerations

When pooling Quantum and 10x libraries on the same flow cell, use a sequencing kit with sufficient read length to accommodate the requirements of both library types. Quantum requires extended Index 1 reads and specific pipeline handling, while 10x requires shorter Index 1 and Read 1 for compatibility with Cell Ranger.

Requirement	Quantum	10x	Both
Read 1	≤82 cycles	28 cycles	82 cycles
Read 2	16 cycles	90 cycles	90 cycles
Index 1	32 cycles	10 cycles	32 cycles
Index 2	8 cycles	10 cycles	10 cycles
			214 cycles

Important:

- Although a 100-cycle kit may be used for Quantum libraries alone, in order to support both Quantum and 10x requirements, a kit with more cycles must be used, as long as Index 1 is at least 32 bases, Index 2 is at least 8 bases, and Read 2 is at least 16 bases.

Generating Separate FASTQ Files


When pooling Quantum and 10x libraries:

- Run FASTQ generation (demultiplexing) twice:
 - Once for Quantum with OverrideCycles:
 - Y82;I8U24;I8N2;Y90
 - Once for 10x with suggested OverrideCycles:
 - Y28N54;I10N22;I10;Y90
- Alternatively, if supported, use BCL Convert with library-specific OverrideCycles settings.
- If possible, sequence libraries on separate flow cell lanes.

Additional tips:

- Work closely with your sequencing core or facility.
- Coordinate with your bioinformatics team to ensure proper demultiplexing and pipeline compatibility.
- Be prepared to handle the demultiplexing and analysis as two distinct workflows.

Figure 1: Parallel library preparation, sequencing, and demultiplexing workflow for QuantumScale and 10x Genomics libraries



Document Revision History

Revision	Revision Date	Document ID	Changes
A	Jun 2025	1352112	Initial release.